


SMITH PETERSON MALLETS



SMITH PETERSON MALLETS

Small	761-9017-0
Large	761-9028-0

TORONTO MALLETS



TORONTO MALLETS

Metal 32 oz.	760-9032-0
Plastic 24 oz.	760-9024-0

MALLET



MALLET

Nylon one side
30mm diameter 10" overall length

761-1030-0

BONE TAMPS



763-9505-0



763-9513-0



763-9510-0

BONE TAMPS 9.5" length 240mm

Rectangular Tip	5mm	763-9505-0
Round Bone Tamp	13mm	763-9513-0
Offset Bone Tamp	10mm	763-9510-0

Offset Bone Tamp excellent for removing osteophytes.

ROUND BONE TAMPS



ROUND BONE TAMPS 9.5" length 240mm

2mm	763-9502-0
3mm	763-9503-0
4mm	763-9504-0
5mm	763-9515-0
6mm	763-9506-0
8mm	763-9508-0



BONE IMPACTORS



BONE IMPACTORS 12" length 310mm

12mm x 6mm 763-1212-0
12mm x 6mm 763-1212-1 SMITH ROBINSON



OFFSET ANGLED BONE TAMP 13" length 330mm

16mm x 5mm 763-1316-1



TAMPS 13.5" length 340mm

12mm x 19mm 763-1312-0

Custom sizes available upon request

LONG IMPACTOR



LONG IMPACTOR
18" length 460mm

10mm round 763-1810-0