

DISK MARKERS



CERVICAL DISK MARKER

The markers easily identify the disk level on the x-ray.
Markers are 10cm in length with a 7mm protective ballstop.
Designed by Edward Dunn, M.D.

Single	790-1070-0
Set of Three	790-1070-1




LUMBAR DISK MARKER

The markers easily identify the disk level on the x-ray.
Markers are 17.5cm in length with a 12.5mm protective ballstop.

Single	790-1013-0
Set of Three	790-1013-1

Saves Time • Reusable • Autoclavable

DISK SPACER SET



DISK SPACER SET

SIZE

9mm	701-1009-0
10mm	701-1010-0
11mm	701-1011-0
12mm	701-1012-0
13mm	701-1013-0
14mm	701-1014-0
16mm	701-1016-0
18mm	701-1018-0
T-Inserter	701-1000-0

"WHITE" RUSSIAN FORCEP



Tamp End

"WHITE" RUSSIAN FORCEP

New forcep designed to facilitate collection, delivery and packing of bone graft. Very useful for bone grafting adjacent to instrumentation, under rods and crosslinks.

Designed by Andrew White, M.D.

8"	780-8005-0
10"	780-1005-0

JEFFERSON FORCEP



JEFFERSON FORCEP

7" length 175mm

1 x 2 Teeth	780-7012-0
Serrated	780-7013-0

LIFE INSTRUMENT ARM WITH UNIVERSAL ATTACHMENT



A

LIFE INSTRUMENT ARM

400mm column 703-1600-0 (includes table clamp)

*Two styles of attachments available.
Must be permanently fixed to arm upon order.*




UNIVERSAL ATTACHMENT (A)

703-1201-0 (shown at left)

LOW PROFILE ATTACHMENT (B)

**Holds Transverse Cervical Retractor Frame
Eliminates Cervical Retraction Movement During Operation**


703-1203-0
(Transverse Retractor Frame not included)



702-0006-D Insert Only

See page 10.4 for outer case and lid information.

LIFE INSTRUMENT L-SHAPE COLUMN WITH ARTICULATING ARM



LIFE INSTRUMENT L-SHAPE COLUMN WITH ARTICULATING ARM

L-Shape Column (450x100mm) with Articulating Arm (400mm extension)

703-1800-L (includes table clamp)

Universal Attachment (sold separately)

703-1201-0

DURA & VESSEL REPAIR INSTRUMENTS



TITANIUM MICRO NEEDLEHOLDER
 12" length 300mm
 Curved with Lock Round Handle for Better Control

782-1201-1



BAYONET MICRO NEEDLEHOLDER
 12" length 300mm
 Curved with Lock

782-1201-3



MICRO SCISSORS
 Stainless Steel
 12" length 300mm
 Curved Round Handle for Better Control

783-1200-1



ATRAUMATIC MICRO FORCEPS
 12" length 300mm
 Debaquey Style 1.2mm Round Handle for Better Control

780-1201-0

BAYONET BIPOLAR FORCEP**BAYONET BIPOLAR FORCEP**
6" working length 1.0mm tip

780-1010-3 Straight Tip

780-1011-3 Angled Tip

KNIFE BLADE HOLDER**KNIFE BLADE HOLDER**
6" working length

705-1010-3

MICRO SCISSOR**MICRO SCISSOR**
7" length

801-0700-0

BONE MILL**BONE MILL**

Includes Cutting Cylinder and Wrench

787-1001-0

Wrench can be sold separately

787-1001-2



TO ORDER, CALL 800-925-2995

JET WIRE TWISTER II



JET WIRE TWISTER II
Stainless Steel

781-1100-1

QUEEN SQUARE PERCUSSION HAMMER



QUEEN SQUARE PERCUSSION HAMMER

Nylon hammer 704-1200-0

LIFE ILLUMINATION SYSTEM

Any Retractor can be modified with an attachment to hold the light bundle.



SOLARC LAMP
24W Light

703-0000-0

CABLE

SIZE

4.5mm	6' length	703-6045-0
5mm	6' length	703-6005-0

LIGHT BUNDLE

SIZE

3mm	180mm	703-0003-0
3mm	150mm	703-0013-0

INSTRUMENTS WITH LONGER WORKING LENGTH

AVAILABLE UPON REQUEST

All instruments can be made with any size handle or stainless steel working end.



DOWN PUSHING CURETTE

15" overall length 7" handle 380mm 180mm handle

719-1510-0S



O'BRIEN CURETTES

15" overall length 7" handle 380mm 180mm handle

8mm 719-1508-0S

11mm 719-1511-0S



CONE RING CURETTES

15" overall length 7" handle 380mm 180mm handle

SIZE	STRAIGHT	ANGLED	LATERAL ANGLED
1	711-1501-0S	711-1501-1S	711-1501-3S
2	711-1502-0S	711-1502-1S	711-1502-3S
3	711-1503-0S	711-1503-1S	711-1503-3S
4	711-1504-0S	711-1504-1S	711-1504-3S

INSTRUMENTS WITH LONGER WORKING LENGTH


AVAILABLE UPON REQUEST



LONG HANDLE CURETTES

15" overall length 7 handle 380mm 180mm handle

SIZE	STRAIGHT	ANGLED	REV. ANGLED
3/0	710-1530-0S	710-1530-1S	710-1530-2S
2/0	710-1520-0S	710-1520-1S	710-1520-2S
0	710-1500-0S	710-1500-1S	710-1500-2S
1	710-1501-0S	710-1501-1S	710-1501-2S
2	710-1502-0S	710-1502-1S	710-1502-2S
3	710-1503-0S	710-1503-1S	710-1503-2S
4	710-1504-0S	710-1504-1S	710-1504-2S
5	710-1505-0S	710-1505-1S	710-1505-2S
6	710-1506-0S	710-1506-1S	710-1506-2S
7	710-1507-0S	710-1507-1S	710-1507-2S



PATTIE TRAY

701-0507-0

CUSTOMS


9" 230mm



All instruments can be customized to surgeons individual preference of handle size and working length.

Handles shown at actual size.

7" 175mm



5" 125mm




5" 125mm



TO ORDER, CALL 800-925-2995

PERSONALIZED COBB ELEVATOR

Great for Gifts!



**All instruments can be personalized on handle
with surgeon's name.**

Gift box included.